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**LAMPIRAN**

**NOMOR IZIN EDAR  
ALAT KESEHATAN  
KEMENKES RI AKL 20304127124**

No.	Deskripsi	Tipe / Kode
1	Power adaptor MIC-PA	
2	Tubes and caps MIC-TUBES	
3	Capping tool MIC-CT	
4	Myra loading block for Mic racks MYRA-LBMIC	

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## micIVD Platform Specifications

### Hardware

Thermal Performance	
Temperature Accuracy	±0.25°C
Temperature Uniformity	±0.1°C
Temperature Input Range	40 – 99°C

Optical Detection System		
Detectors	Photodiode	
Sources	LED	
Green Channel	Ex 465nm	Em 510nm
Yellow Channel	Ex 540nm	Em 570nm
Orange Channel	Ex 585nm	Em 618nm
Red Channel	Ex 635nm	Em 675nm

Reaction Tubes	
Samples per instrument	48
Reaction Volume Range	5 – 30µL

Operating Environment	
Ambient Temperature	18 – 30°C
Relative Humidity	20 – 80%

Physical	
Dimensions	150x150x130mm
Cover Panel Colour	BMS Red
Lid Insert Panel	Mic IVD

## Software

micIVD software to be based on micPCR v2.10 public release software with the following changes:

Device Identification	
Unique Device ID	micIVD software is locked to work with micIVD device ID instruments only. micPCR RUO micPCR software will not detect micIVD devices

Branding	
Splash screen	micIVD specific
Software title	micIVD vX.Y.Z (<User level>) where X.Y.Z is the version number and <User level> is the user permission level
About screen	micIVD specific
User manual	micIVD specific

Access Control	
Levels provided	4 levels <ul style="list-style-type: none"> <li>• Technician</li> <li>• Analyst</li> <li>• Administrator</li> <li>• Maintenance</li> </ul>
Maintenance	Only accessible for IVD distributors. Full control of software and in addition can <ul style="list-style-type: none"> <li>• Can perform TVS temperature verification</li> <li>• Select by using special 'micIVD Maintenance' shortcut created by 'micIVD Maintenance' installer</li> </ul>
Administrator	Full software control of software including <ul style="list-style-type: none"> <li>• Creating and editing of assays and templates</li> <li>• Editing of run profiles</li> <li>• Editing of instrument properties and updating firmware</li> <li>• Cannot perform TVS temperature verification</li> <li>• Default user level if no micIVD Windows user groups defined</li> </ul>
Analyst	All Technician permissions, and in addition <ul style="list-style-type: none"> <li>• Can create and edit analyses and reports</li> <li>• Can view run and assay profiles</li> </ul>
Technician	<ul style="list-style-type: none"> <li>• Create runs from templates</li> <li>• Can edit sample information</li> <li>• Can start and abort runs</li> <li>• Can only view (not create or edit) analyses or reports</li> <li>• Cannot view run or assay profiles</li> </ul>

Analysis Types	
Included	<ul style="list-style-type: none"> <li>• Cycling</li> </ul>

	<ul style="list-style-type: none"> <li>• Melt</li> <li>• Absolute Quantification</li> </ul>
Excluded	<ul style="list-style-type: none"> <li>• High Resolution Melt</li> <li>• Identifier</li> <li>• Allelic Discrimination</li> <li>• Relative Quantification</li> <li>• Projects</li> </ul>

<b>File Handling</b>	
Identity and protection	<ul style="list-style-type: none"> <li>• All files include Unique ID to lock file to micIVD software</li> <li>• All files include digital signature to prevent unauthorised modifications</li> </ul>
Assay File	<ul style="list-style-type: none"> <li>• Created and edited in Administrator Mode</li> <li>• Includes thermal profile and reaction and analysis information. Added to runs to shared information.</li> <li>• Technician and Analyst Mode users can only add assays with matching thermal profiles to runs</li> <li>• In Administrator and Maintenance Mode assay thermal profile will be merged with existing run thermal profile when assay is added</li> </ul>
Run Template File	<ul style="list-style-type: none"> <li>• Created/controlled by Administrator Mode users for use by Technician and Analyst users to create runs</li> <li>• Installed to a fixed directory (refer 'Install Template' menu below)</li> <li>• Can include additional predefined data (e.g. controls, sample names, run information)</li> <li>• New thermal profile, or updates to existing thermal profile requires a new run template file to be installed by user (refer 'Install Template' menu below).</li> </ul>
Install Template menu	In Technician and Analyst Modes a new menu item 'Install Template' will be added to the 'Help' menu to allow users to install updated templates provided by Administrator.
New Run	In Technician and Analyst Modes will automatically use installed template file. Prompt for selection if more than one template installed. In Administrator and Maintenance Modes will create an empty run.
New Run From Template	Hidden in Technician and Analyst Modes
New Assay	Hidden in Technician and Analyst Modes
New Project	Removed
Save as Template	Hidden in Technician and Analyst Modes
Template Filename Extension	micivdtemplate
Run Filename Extension	micivdrun
Assay Filename Extension	micivdassay

<b>Assays</b>	
Automatic Analysis Generation (per target)	Setting to automatically add analysis at conclusion of successful run. All data must be present (not aborted or failed) and not suspected of optical failure

<b>Sample Setup</b>	
Available assays	Assay list on sample editor shows all assays included in template file.
Sample Type Options	Unknown, Standard, PositiveControl, NegativeControl, NTC, NRT
Groups	Groups columns and setup area removed (only required for Relative Quantitation Analysis type)
Available columns	Colour, Name, Type, Assay, Standards Concentration, Input DNA Concentration, Input RNA Concentration, RIN, Comment
Removed columns	Groups, Multiplex Standards, Reference Material
Editing period	Sample setup cannot be changed after run starts in Technician or Analyst Mode

<b>Analysis and Reporting</b>	
Automatic report generation	A report will be automatically at the end of a completed run. The analysis report for each target will be listed in the report configuration panel. All standard report properties will be included.
Analysis parameters settings	Disabled (non editable) in Technician Mode
Analysis parameters save/load to assay	Hidden in Technician Mode Save to assay not hidden in Analysis Mode
Sample Selector	In Technician Mode cannot change sample selection but can hide/show samples. In Analyst, Administrator and Maintenance Mode can enable/disable samples, and hide/show samples.
Sample Selector Grouping	Remove group by 'Sample Group'

<b>Temperature Verification System (TVS)</b>	
Verification of Mic Temperature Sensor	Wizard driven system to verify Mic temperature accuracy.
Traceability	TVS calibration traceable to NATA reference standard.
Availability	Only available in Maintenance Mode

## MIC IVD Labels Artwork

Mic IVD label is placed in the position on the carton as shown below:

**MIC IVD Carton (39x39x26cm) (weighs 6kg).**